

# NOVAMED<sup>®</sup>

*Your Health is important to us*



*Orthopedics  
Product Catalog*



[www.novamed.tc](http://www.novamed.tc)

# About Us

*A leading medical company in Turkish market; NOVAMED founded in 2004 in medical sector. We manufacture medical products and marketing .*

*We focused especially;*

- **Orthopedic Products***
- **Medical Stockings***
- **Silicone Products***
- **Visco Pillows***

*Our company headquarter located in Istanbul, Turkey. Our product range is covering all doctors, hospitals, surgery centers, pharmacy shops, physician's offices, sub dealers and patients' needs with assurance of good quality products and good quality service conception.*

*NOVAMED has CE certificate with proper 93/42/EEC European Community Health Law conditions, and ISO 13485:2012 certificate. Also we have registered trade mark "NOVAMED".*

*Our company participate in major medical fairs. As NOVAMED, our products are exported more than 30 countries.*

*Your Health  
is important  
to us*

*Index*

---

# alth tant

---

<i>Material Preferences</i>	<b>4</b>	<i>Calf and Shorts</i>	<b>63-65</b>
<i>Quality Policy &amp; Certificates</i>	<b>5</b>	<i>Elastic Products</i>	<b>66</b>
<i>Collar</i>	<b>7-10</b>	<i>Pillows</i>	<b>67-69</b>
<i>Corset</i>	<b>11-22</b>	<i>Medical Stockings</i>	<b>70-75</b>
<i>Shoulder and Arm</i>	<b>23-27</b>	<i>Lymph Edema Arm Sleeve</i>	<b>76-78</b>
<i>Elbow</i>	<b>28-29</b>	<i>Silicone Products</i>	<b>79-91</b>
<i>Hand and Wrist</i>	<b>30-41</b>	<i>Support Products</i>	<b>92-96</b>
<i>Knee</i>	<b>42-51</b>	<i>Other Products</i>	<b>97-98</b>
<i>Ankle and Foot</i>	<b>52-62</b>		

# The Exclusive Advantages Of The Novamed Products




Neoprene; It is laminated with jersey and elastic terry cloth. Neoprene material has high durability and also provides ease of movement, and comfort heat by its elastic structure. Increased body temperature accelerates recovery process by improving the blood circulation.

Between velcro sticky fabric and sweat absorbent cotton surface is sponge. Sponge fabric with high density provides the require stabilization by wrapping the body



It's flexible, and high air permeable cotton fabric. Flexible design wrap body to provide required support and comfort.



It's flexible, high air permable, moisture proof, and anatomically shaped special knitted fabric. With special knitting design it provives support for requirement parts of the body and does not the generate the pressure and tightnes feeling over joints zone. Fits perfectly to body and provides maximum comfort by its massage effects.







**ORTHOPEDIC PRODUCTS**

*Novamed Orthopedic Products used for many purposes such as controlling an organ or joint in a certain part of the body, restricting movement in the specified direction, helping the movement in general, distributing the body burden evenly, helping rehabilitation in the recovery process of fractures or reducing pain.*



**CODE : 1700**

**CERVICAL COLLAR WITH CHIN SUPPORT**

Used for cervical degenerations, traumas, after operation at neck arthritis and neck hernia. Produce from plastic material. It has chin support. Height and length are adjustable with velcros. Designed according to anatomical shape of neck.



Neck Circ. (cm)	22-27	27-32	32-38	38-44	44-50
Size	S	M	L	XL	XXL
Article Nr.	1700-1	1700-2	1700-3	1700-4	1700-5



**CODE : 1800**

**CERVICAL COLLAR WITHOUT CHIN SUPPORT**

Used for cervical degenerations, traumas, after operations. Produce from plastic material. Height and length are adjustable with velcros. Designed according to anatomical shape of neck.



Neck Circ. (cm)	22-27	27-32	32-38	38-44	44-50
Size	S	M	L	XL	XXL
Article Nr.	1800-1	1800-2	1800-3	1800-4	1800-5



**CODE : 1300**

**PHILADELPHIA CERVICAL COLLAR**

Used for neck fractures, cervical degenerations and traumas, after operation and in urgent patient transportation. Produced as a two part and made of antiallergic plastic. Designed according to anatomic shape of neck.



Neck Circ. (cm)	25-31	32-38	39-46
Size	S	M	L
Article Nr.	1300-S	1300-M	1300-L



**CODE : 1350**

**PHILADELPHIA CERVICAL COLLAR FOR TRACHEOSTOMY**

Used for tracheostomy patients, cervical degenerations and traumas also for patient transport and after operation. Produced as a two part and made of antiallergic plastic. Trace opening enables to treat neck area in urgent cases. Designed according to anatomic shape of neck.



Neck Circ. (cm)	25-31	32-38	39-46
Size	S	M	L
Article Nr.	1350-S	1350-M	1350-L



**CODE : 2400**

**SPONGE CERVICAL COLLAR**

*Used for cervical degenerations, trauma, calcification, mild pain and after operation.  
Designed according to anatomic shape of neck.*



Neck Circ. (cm)	22-27	27-32	32-38	38-44	44-50
Size	S	M	L	XL	XXL
Article Nr.	2400-1	2400-2	2400-3	2400-4	2400-5



**CODE : 2410**

**SPONGE CERVICAL COLLAR**

*Used for cervical degenerations, trauma, calcification, mild pain and after operation.  
Designed according to anatomic shape of neck.*



Neck Circ. (cm)	22-27	27-32	32-38	38-44	44-50
Size	S	M	L	XL	XXL
Article Nr.	2410-1	2410-2	2410-3	2410-4	2410-5



**CODE : 1500**

**NELSON CERVICAL COLLAR**

Used for cervical degenerations, traumas, and soft tissue damages, transportation of patients, and after operations. Designed according to anatomic shape of neck.



Neck Circ. (cm)	25-28	29-33	38-42	39-46
Size	0	1	2	3
Article Nr.	1500-0	1500-1	1500-2	1500-3



**CODE : 1325**

**WIRE FRAME CERVICAL COLLAR**

Used for cervical degenerations, severe traumas, after operation, at neck arthritis and neck hernia. Designed according to anatomic shape of neck. Height and length are adjustable.



Size	Small	Regular
Article Nr.	1325-S	1325-L



**CODE : 4601**

**KNITTED LUMBOSACRAL CORSET**

Used for lumbosacral traumas, lomberdisc, hernias, lomber spondylosis, and after operation. Consists four flexible steel on back side which provides perfect support to lumbar to relieves the pain in muscles. Produced from breathable, flexible knitted fabric.



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4601-S	4601-M	4601-L	4601-XL	4601-XXL



**CODE : 4602**

**KNITTED LUMBOSACRAL CORSET WITH EXTRA SUPPORT**

Used for lumbosacral traumas, lomber disc, hernias, lomber spondylosis, and after operations. Produced from breathable, flexible knitted fabric. Four flexible steel stays fit according to shape of body back. Extra cross straps offers additional support and stabilisation to lumbar area.



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4602-S	4602-M	4602-L	4602-XL	4602-XXL




**CODE : 4603**

**KNITTED LUMBOSACRAL CORSET WITH SUPPORT**

Used for lumbosacral traumas, lumbar disc, hernias, lumbar spondylosis, and after operations.. Produced from breathable, flexible knitted fabric. Four flexible steel stays fit according to shape of body back. Suitable for underneath of daily clothing.



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4603-S	4603-M	4603-L	4603-XL	4603-XXL



**CODE : 4604**

**KNITTED LUMBOSACRAL CORSET WITH SILICONE SUPPORT**

Used for lumbosacral traumas, lumbar disc, hernias, lumbar spondylosis, and after operation. Consists four flexible steel on back side which provides perfect support to lumbar to relieve the pain in muscles. Accelerates blood circulation, and offers massage effect with silicone pad in lumbar area that relieves the pain. Produced from breathable, flexible knitted fabric.




Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4604-S	4604-M	4604-L	4604-XL	4604-XXL



**CODE : 7700**

**HYPEREXTENSION CORSET**

*Produced from special aluminium. Parts where touches to skin is produced from terry cloth. Restricts the flexion of the body and offers hyperextension position. Restricts spinal movement in certain point. Provides anatomical posture and aids recovery process.*



<b>Body Circ. (cm)</b>	<b>40-49</b>	<b>50-60</b>
<b>Size</b>	<b>S</b>	<b>L</b>
<b>Article Nr.</b>	<b>7700-S</b>	<b>7700-L</b>



**CODE : 7710**

**HYPEREXTENSION CORSET**

*Produced from special aluminium. Parts where touches to skin is covered with leather. Length and width adjustment is possible by the screw system on the orthosis. Restricts the flexion of the body and offers hyperextension position. Restricts spinal movement in certain point. Provides anatomical posture and aids recovery process.*



<b>Body Circ. (cm)</b>	<b>40-46</b>	<b>47-53</b>	<b>54-60</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>7710-S</b>	<b>7710-M</b>	<b>7710-L</b>



**CODE : 5600**

**DORSOLUMBAR CORSET 26 CM**

Used for compression fractures, osteoporosis, arthritis, lordosis, scoliosis and after operations. Provides maximum spinal coverage in required situations. Produce from breathable and sweat absorbing corset fabric. There are two pieces of plastic balen in both right and left sides. Also there are two steel stays on back side which provide perfect support. Steel stays can be removable if needed.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	5600-S	5600-M	5600-L	5600-XL	5600-XXL



**CODE : 5610**

**DORSOLUMBAR CORSET 32 CM**

Used for compression fractures, osteoporosis, arthritis, lordosis, scoliosis and after operations. Provides maximum spinal coverage in required situations. Produce from breathable and sweat absorbing corset fabric. There are two pieces of plastic balen in both right and left sides. Also there are two steel stays on back side which provide perfect support. Steel stays can be removable if needed.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	5610-S	5610-M	5610-L	5610-XL	5610-XXL



**CODE : 2500**

**POSTUREX CORSET**

Used for correction of posture disorders and keeps the body upright. Used in young people, and employees who suffers postural disorders who spent a lot of time in desk job. Produced from breathable, and flexible fabric. Used in young people, and employees who suffers postural disorders who spent a lot of time in desk job. Produced from breathable, and flexible corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	2500-S	2500-M	2500-L	2500-XL	2500-XXL



**CODE : 5700**

**ABDOMINAL CORSET**

Used after operation, and preventing abdomen prolapse after birth. Produced from breathable, and flexible corset fabric.

Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	5700-S	5700-M	5700-L	5700-XL	5700-XXL



**CODE : 4600**

**LUMBOSACRAL CORSET 26 CM**

Used for lumbosacral traumas, lomber disc, hernias, lomber spondylodiscitis, and after lumbar operations. Plastic stays on back side, and steel stays on left and right side of spine which removable. Provide extra stabilisation to lumbar area with cross straps. Produced from breathable, and flexible corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4600-S	4600-M	4600-L	4600-XL	4600-XXL



**CODE : 4610**

**LUMBOSACRAL CORSET 32 CM**

Used for lumbosacral traumas, lomber disc, hernias, lomber spondylodiscitis, and after lumbar operations. Plastic stays on back side, and steel stays on left and right side of spine which removable. Provide extra stabilisation to lumbar area with cross straps. Produced from breathable, and flexible 32cm corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4610-S	4610-M	4610-L	4610-XL	4610-XXL



**CODE : 5400**

**LUMBOSTAD CORSET**

Used for chronic lumbago, sciatica, waist hernia, after given birth, support for back and abdominal at operations and after other surgical operations. Four plastic stays on back side provide perfect support for spine. Produced from breathable, and flexible corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	5400-S	5400-M	5400-L	5400-XL	5400-XXL



**CODE : 5800**

**COLOSTOMY CORSET**

Produced from breathable, flexible corset fabric. Suitable for daily activities. Support for abdominal area. Decreases the risk of re-occurring hernia after operations. Suitable for underneath clothing.

Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	5800-S	5800-M	5800-L	5800-XL	5800-XXL



**CODE : 6600**

**HERNIA CORSET**

Produced from breathable, flexible, and sweat absorbing corset fabric. Offers anatomically designed pad which applies compression to herniated area. Usage of the umbilical hernia corset is not interfere with daily activities, and it is suitable for underneath wearing in daily clothing.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	6600-S	6600-M	6600-L	6600-XL	6600-XXL



**CODE : 5701**

**PREGNANCY CORSET**


Produced from breathable, flexible, and sweat absorbing corset fabric. Designed for before and after giving birth. Provides support in pregnancy period, and prevents lumbar pain and back aches. Prevents sagging after birth. Suitable for daily usage, and activities.



Colors



Waist Circ. (cm)	85-95	95-105	105-115	115-125	125-135
Size	S	M	L	XL	XXL
Article Nr.	5701-S	5701-M	5701-L	5701-XL	5701-XXL




**CODE : 4605**

**LUMBAR BACK WITH DOUBLE SUPPORT**

*Produced from neoprene fabric. double support on backside provides good support to both sides of lumbar area. Offers individualized, adjustable pressure levels by straps. Provides proper heating during activity.*



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4605-S	4605-M	4605-L	4605-XL	4605-XXL



**CODE : 4606**

**LUMBAR BACK SUPPORT WITH PAD**

*Design is made from neoprene fabric, Backside soft pad provides perfect support. Design offers individualized, adjustable pressure levels by straps. Fabrication provides proper heating during activity. It should be weared around waist.*



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	4606-S	4606-M	4606-L	4606-XL	4606-XXL



**CODE : 1100**

**WOOL WAIST CORSET**

*Used for rheumatism, waist pains, keeping abdomen muscle warm. Produced from elastic wool fabric.*



<b>Waist Circ. (cm)</b>	<b>70-80</b>	<b>80-90</b>	<b>90-100</b>	<b>100-110</b>	<b>110-120</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1100-S</b>	<b>1100-M</b>	<b>1100-L</b>	<b>1100-XL</b>	<b>1100-XXL</b>



**CODE : 1160**

**WOOL WAIST CORSET WITH VELCRO**

*Used for rheumatism, waist pains, keeping abdomen muscle warm. Provides comfort usage with velcros. Produced from elastic wool fabric.*



<b>Waist Circ. (cm)</b>	<b>70-80</b>	<b>80-90</b>	<b>90-100</b>	<b>100-110</b>	<b>110-120</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1160-S</b>	<b>1160-M</b>	<b>1160-L</b>	<b>1160-XL</b>	<b>1160-XXL</b>



**CODE : 3900**

**CHEST CORSET**

Used for after heart operation. Produced from breathable, and sweat absorbing material.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	3900-S	3900-M	3900-L	3900-XL	3900-XXL



**CODE : 3925**

**CARDIOLOGY CORSET**

Used for after heart operation. Fits around the upper rib cage. Provides compression and immobilization for the upper chest area. Ideal for support and minimizing movement following rib and sternum fractures. Produced from breathable, and sweat absorbing corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	3925-S	3925-M	3925-L	3925-XL	3925-XXL



**CODE : 3930**

**CHEST CORSET**


Used for after heart operation. Fits around the upper rib cage. Provides compression and immobilization for the upper chest area. Ideal for support and minimizing movement following rib and sternum fractures. Produced from breathable, and sweat absorbing corset fabric.



Colors



Waist Circ. (cm)	70-80	80-90	90-100	100-110	110-120
Size	S	M	L	XL	XXL
Article Nr.	3930-S	3930-M	3930-L	3930-XL	3930-XXL



**CODE : 4100**

**HERNIA BELT SINGLE**

*Designed to push back to the herniated tissue. Suitable for daily use. Prevent reoccurs of hernia after operations, and provides support to healing process. Produced from flexible cotton fabric.*



<b>Size</b>	<b>Right</b>	<b>Left / Standard</b>
<b>Article Nr.</b>	<b>4100-R</b>	<b>4100-L</b>




**CODE : 4200**

**HERNIA BELT DOUBLE**

*Designed to push back to the herniated tissue. Suitable for daily use. Prevent reoccurs of hernia after operations, and provides support to healing process. Produced from flexible cotton fabric.*



<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>4200</b>




**CODE : 300**

**CLAVICLE BANDAGE**

*Used for shoulder fractures, stabilisation of shoulder, correction of posture disorder, and keeping the body up. Produced from sponge fabric.*



<b>Arm Circ. (cm)</b>	<b>26-31</b>	<b>32-37</b>	<b>38-44</b>	<b>45-51</b>	<b>52-58</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>300-S</b>	<b>300-M</b>	<b>300-L</b>	<b>300-XL</b>	<b>300-XXL</b>




**CODE : 305**

**CLAVICLE BANDAGE ADJUSTABLE**

*Used for shoulder fractures, stabilisation of shoulder, correction of posture disorder, and keeping the body up. Provides correct posture, and needed tension for body with its adjustable straps. Produced from sponge fabric.*



<b>Arm Circ. (cm)</b>	<b>26-31</b>	<b>32-37</b>	<b>38-44</b>	<b>45-51</b>	<b>52-58</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>305-S</b>	<b>305-M</b>	<b>305-L</b>	<b>305-XL</b>	<b>305-XXL</b>



**CODE : 7800**

**SHOULDER PAD**

Used for after shoulder traumas, before and after surgery, shoulder bursitis, and shoulder injuries. Produced from antiallergic and antibacterial neoprene fabric.



Biceps Circ. (cm)	22-24	24-27	27-30	30-34	34-38
Size	S	M	L	XL	XXL
Article Nr.	7800-S	7800-M	7800-L	7800-XL	7800-XXL



**CODE : NV.7801**

**KNITTED SHOULDER PAD**

Used for after shoulder traumas, before and after surgery, shoulder bursitis, and shoulder injuries. Produced from breathable, flexible, antiallergic and antibacterial knitted fabric.



Biceps Circ. (cm)	22-24	24-27	27-30	30-34	34-38
Size	S	M	L	XL	XXL
Article Nr.	NV.7801-S	NV.7801-M	NV.7801-L	NV.7801-XL	NV.7801-XXL



**CODE : 1450**

**ARM SLING WITH ADJUSTABLE ANGLE**

*Produced from sponge fabric. Used for soft tissue damages. Provides stabilisation when arm requires to remain at 30°. Offers support for recovery process.*




Arm Circ. (cm)	22-25	25-28	28-31	31-34	34-37
Size	S	M	L	XL	XXL
Article Nr.	1450-S	1450-M	1450-L	1450-XL	1450-XXL




**CODE : 3800**

**VELPEAU BANDAGE**

*Produced from sponge fabric. Used for humeral fractures, shoulder luxation, and after operations. Offers stabilisation with wraparound belt and provides support for healing process.*



Arm Circ. (cm)	22-25	25-28	28-31	31-34	34-37
Size	S	M	L	XL	XXL
Article Nr.	3800-S	3800-M	3800-L	3800-XL	3800-XXL



**CODE : 3200**

**ARM SLING WITH SPONGE**

Used for soft tissue damages that require the arm to remain motionless. Used after operation for arm sling when necessary. Produced sponge fabric.




Arm Circ. (cm)	22-25	25-28	28-31	31-34	34-37
Size	S	M	L	XL	XXL
Article Nr.	3200-S	3200-M	3200-L	3200-XL	3200-XXL




**CODE : 3250**

**ARM SLING WITH MESH**

Used for soft tissue damages that require the arm to remain motionless. Used after operation for arm sling when necessary. Produced mesh fabric.



Arm Circ. (cm)	22-25	25-28	28-31	31-34	34-37
Size	S	M	L	XL	XXL
Article Nr.	3250-1	3250-2	3250-3	3250-4	3250-5



**CODE : 3260**

**ARM SLING CLASSIC**

*Produced from velour fabric. Offers individual adjustment by cutting on needed lenght. Available for both right and left side. Presents friendly-use to patient with natural material.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>3260</b>



**CODE : 6720**

**EPICONDYLITIS BANDAGE WITH SUPPORT**

*Used for elbow traumas, elbow inflammation, ligament injuries, healing process after surgery, and protect in sports activities. Offers support with special designed silicone pad and relieves the pain in tendon muscles. Produced from antiallergic and antibacterial neoprene fabric.*



<b>Elbow Circ. (cm)</b>	<b>22-24</b>	<b>24-26</b>	<b>26-28</b>	<b>28-30</b>	<b>30-32</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6720-S</b>	<b>6720-M</b>	<b>6720-L</b>	<b>6720-XL</b>	<b>6720-XXL</b>



**CODE : 6710**

**EPICONDYLITIS BANDAGE (LONG)**

Used for elbow traumas, elbow inflammation, ligament injuries, healing process after surgery, and protect in sports activities. Offers extra support with special designed pad and straps. Produced from antiallergic and antibacterial neoprene fabric.



Elbow Circ. (cm)	22-24	24-26	26-28	28-30	30-32
Size	S	M	L	XL	XXL
Article Nr.	6710-S	6710-M	6710-L	6710-XL	6710-XXL



**CODE : 6700**

**EPICONDYLITIS BANDAGE (SHORT)**

Used for elbow traumas, elbow inflammation, ligament injuries, healing process after surgery, and protect in sports activities. Offers support with special designed pad and relieves the pain in tendon muscles. Produced from antiallergic and antibacterial neoprene fabric.


Size	Standard
Article Nr.	6700



**CODE : 6725**

**ARM ORTHOSIS WITH ADJUSTABLE ANGLE**

Produced from terry cloth where the parts touches the skin Surface covered with plastic material. Angle adjustment is possible with flawless mechanism. Additional supports with straps from shoulder. Options for adjustment for length and height.



Arm Circ. (cm)	30-37	38-45		
Size	S	L		
Article Nr.	6725-1	6725-3	6725-S1	6725-S3



**CODE : NV.7301**

**KNITTED ELBOW**

Designed to improve blood circulation, and accelerate recovery process. Provides compression on muscles with micromassage effect pads on both sides of elbow. Comfortable for using in all sport activities as design avoids irritation. Produced from breathable, and flexible knitted fabric.



Elbow Circ. (cm)	22-24	24-26	26-28	28-30	30-32
Size	S	M	L	XL	XXL
Article Nr.	NV.7301-S	NV.7301-M	NV.7301-L	NV.7301-XL	NV.7301-XXL



**CODE : 1000**

**HAND WRIST SPLINT**

Used for hand and ankle traumas, sprains, and after operation to keep the hand in rest position. Support for stabilisation of wrist. Produced from breathable, and flexible fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	1000-1 1000-S1	1000-2 1000-S2	1000-3 1000-S3	1000-4 1000-S4	1000-5 1000-S5



**CODE : 2800**

**HAND WRIST SPLINT (LONG)**

Used for hand and ankle traumas, sprains, and after operation to keep the hand in rest position. Provides support with aluminium shaped as spoon for stabilisation of wrist. Produced from breathable, and flexible fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	2800-1 2800-S1	2800-2 2800-S2	2800-3 2800-S3	2800-4 2800-S4	2800-5 2800-S5



**CODE : 7200**

**HAND WRIST SPLINT (NEOPRENE)**

Used for hand and ankle traumas, sprains, and after operation to keep the hand in rest position. Support for stabilisation of wrist. Produced from antiallergic, and antibacterial neoprene fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7200-1 7200-S1	7200-2 7200-S2	7200-3 7200-S3	7200-4 7200-S4	7200-5 7200-S5



**CODE : 7250**

**HAND WRIST SPLINT WITH THUMB SUPPORT (NEOPRENE)**

Used for when wrist and thumb require rest. Has adjustable aluminium support for restricting movement of hand and thumb. Produced from antiallergic and antibacterial neoprene fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7250-1 7250-S1	7250-2 7250-S2	7250-3 7250-S3	7250-4 7250-S4	7250-5 7250-S5



**CODE : 1900**

**HAND WRIST SPLINT WITH THUMB SUPPORT**

Used for when wrist and thumbs require rest. Has adjustable aluminium support for restricting movement of hand and thumb. Produced from breathable, and flexible fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	1900-1 1900-S1	1900-2 1900-S2	1900-3 1900-S3	1900-4 1900-S4	1900-5 1900-S5



**CODE : 7180**

**THUMB SUPPORT**

Used for when thumb require rest. Offers stabilisation with two stays for perfect fit. Adjustable straps provides individual length. Produced from breathable, antibacterial, antiallergic material.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7180-1 7180-S1	7180-2 7180-S2	7180-3 7180-S3	7180-4 7180-S4	7180-5 7180-S5



**CODE : 1925**

**RESTING SPLINT**

Produced from thermoplastic material, and terry cloth fabric. Provides support for hand and wrist when needed. Used for stabilisation of hand and wrist after operations.



Wrist Circ. (cm)	14-18	18-22		
Size	S	L		
Article Nr.	1925-1	1925-3	1925-S1	1925-S3



**CODE : 7100**

**BASIC HAND WRIST SPLINT**

Used for wrist traumas, sprains, and wrist pains. Used for protection in sport activities. Produced from antiallergic and antibacterial neoprene fabric.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7100-S	7100-M	7100-L	7100-XL	7100-XXL



**CODE : 1003**

**WRIST SPLINT WITH SUPPORT**

*Provides stabilise to the wrist joint. Relieves pain and correct instability in wrist. Offers movement restriction of wrist and helps healing process. Adjustable straps offers individual fit, and comfort.*



<b>Wrist Circ. (cm)</b>	<b>14-16</b>	<b>16-18</b>	<b>18-20</b>	<b>20-22</b>	<b>22-24</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1003-1</b> <b>1003-S1</b>	<b>1003-2</b> <b>1003-S2</b>	<b>1003-3</b> <b>1003-S3</b>	<b>1003-4</b> <b>1003-S4</b>	<b>1003-5</b> <b>1003-S5</b>



**CODE : 1006**

**HAND WRIST SPLINT WITH SUPPORT**

*Ensures perfect support and protection to the wrist joint without restricting the movement of the hand. Spoon shaped plastic splint and plastic stays on front provides stabilisation. Level of pressure can be optimised by an easy-to-apply elastic band. Offers easy usage.*



<b>Wrist Circ. (cm)</b>	<b>14-16</b>	<b>16-18</b>	<b>18-20</b>	<b>20-22</b>	<b>22-24</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1006-1</b> <b>1006-S1</b>	<b>1006-2</b> <b>1006-S2</b>	<b>1006-3</b> <b>1006-S3</b>	<b>1006-4</b> <b>1006-S4</b>	<b>1006-5</b> <b>1006-S5</b>



**CODE : 1012**

**HAND WRIST SPLINT CLASSIC**

*Produced from breathable, and flexible fabric. Offers support with spoon shaped plastic stays to the wrist. Offers perfect fit with velcro fasteners.*



<b>Wrist Circ. (cm)</b>	<b>14-16</b>	<b>16-18</b>	<b>18-20</b>	<b>20-22</b>	<b>22-24</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1012-1</b> <b>1012-S1</b>	<b>1012-2</b> <b>1012-S2</b>	<b>1012-3</b> <b>1012-S3</b>	<b>1012-4</b> <b>1012-S4</b>	<b>1012-5</b> <b>1012-S5</b>



**CODE : 1013**

**HAND WRIST SPLINT CLASSIC (LONG)**

*Produced from breathable, and flexible fabric. Offers fast healing process by its long design, and support to the wrist. Presents stabilisation for wrist and hand with spoon shaped plastics stays. Offers perfect fit with velcro fastener.*



<b>Wrist Circ. (cm)</b>	<b>14-16</b>	<b>16-18</b>	<b>18-20</b>	<b>20-22</b>	<b>22-24</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>1013-1</b> <b>1013-S1</b>	<b>1013-2</b> <b>1013-S2</b>	<b>1013-3</b> <b>1013-S3</b>	<b>1013-4</b> <b>1013-S4</b>	<b>1013-5</b> <b>1013-S5</b>



**CODE : 1005**

**THUMB SPLINT**

Produced from breathable sponge fabric. Offers individual stabilisation with straps in wrist area. Designed as anatomically shaped of hand, and offers supports with shapeable aluminium splint. Right and left options available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	1005-1 1005-S1	1005-2 1005-S2	1005-3 1005-S3	1005-4 1005-S4	1005-5 1005-S5



**CODE : 1001**

**KNITTED HAND WRIST SPLINT**

Produced from breathable, flexible, and sweat absorbing knitted fabric. Protects and relieves the veins and nerves with silicone support around ulnar steloid. Plastic L shaped support is anatomically formed to fit to the wrist. Bandage is manufactured for maximum fit, comfort, and aids support injuries, increase blood circulation, and keep body heat and accelaretes the healing.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	1001-1 1001-S1	1001-2 1001-S2	1001-3 1001-S3	1001-4 1001-S4	1001-5 1001-S5



**CODE : 9101**

**WRIST BANDAGE STANDARD**

*Designed to increase the blood circulation and avoid mild stretches in the hand or forearm. Used for pain relieving, strains and inflammations. Produced from breathable, and flexible fabric.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>9101</b>



**CODE : 9102**

**HAND BANDAGE CLASSIC**

*Designed to increase the blood circulation and avoid mild stretches in the hand. Used for pain relieving, strains and inflammations. Produced from breathable, and flexible fabric.*



<b>Wrist Circ. (cm)</b>	<b>14-16</b>	<b>16-18</b>	<b>18-20</b>	<b>20-22</b>	<b>22-24</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>9102-1</b>	<b>9102-2</b>	<b>9102-3</b>	<b>9102-4</b>	<b>9102-5</b>
	<b>9102-S1</b>	<b>9102-S2</b>	<b>9102-S3</b>	<b>9102-S4</b>	<b>9102-S5</b>



**CODE : 7150**

**THUMB SPLINT WITH SUPPORT**

Used for thumb, muscle, ligament strain, and when thumb requires rest. Provides support with aluminium stays to thumb and restricts the movement of thumb. Produced from antiallergic, antibacterial neoprene fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7150-1 7150-S1	7150-2 7150-S2	7150-3 7150-S3	7150-4 7150-S4	7150-5 7150-S5



**CODE : 7170**

**THUMB SPLINT WITH SUPPORT ADVANCED**

Used for thumb, muscle, ligament strain, and when thumb requires rest. Provides support with aluminium stays to thumb and restricts the movement of thumb. Produced from antiallergic, antibacterial neoprene fabric. Right and left sizes options are available.



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	7170-1 7170-S1	7170-2 7170-S2	7170-3 7170-S3	7170-4 7170-S4	7170-5 7170-S5



**CODE : 1004**

**THUMB SUPPORT ELEGANCE**

*Offers stabilisation for thumb, and protects thumb joint. Relieves pain in cases of damages. Used for thumb sprains middle thumb joint damages. Produced from antiallergic, antibacterial fabric.*



Wrist Circ. (cm)	14-16	16-18	18-20	20-22	22-24
Size	S	M	L	XL	XXL
Article Nr.	1004-1 1004-S1	1004-2 1004-S2	1004-3 1004-S3	1004-4 1004-S4	1004-5 1004-S5



**CODE : 1015**

**THUMB SUPPORT ADVANCED**

*Offers stabilisation to the thumb joint. Help healing process of thumb joint damages. Reduces pain in thumb area with the optimal support. Produced antiallergic, antibacterial fabric.*



Wrist Circ. (cm)	14-18	18-22	22-26
Size	S	M	L
Article Nr.	1015-1 1015-S1	1015-2 1015-S2	1015-3 1015-S3



**CODE : 150**

**HAND THERAPY BALL**

Used for activating hand muscles and for hand exercise.  
Produced from high quality silicone.

Size	Green (Soft)	Red (Semi rigid)	Blue (Hard)
Article Nr.	150-1	150-2	150-3



**CODE : 6000**

**MALLET FINGER**

Used for finger fractures, deformities and immobilize to finger in sport injuries. There are holes on the splint and nail side is open.



Finger Circ. (cm)	4-4,5	4,5-5	5-5,5	5,5-6	6-6,5	6,5-7	7-7,5
Size	1	2	3	4	5	6	7
Article Nr.	6000-1	6000-2	6000-3	6000-4	6000-5	6000-6	6000-7




**CODE : 3050**

**BASEBALL SPLINT**

*Used for finger fractures and deformities, sports injuries, and detect finger, and when parts of the finger need rest.*



<b>Finger Circ. (cm)</b>	4,5-5	5-6	6-7
<b>Size</b>	S	M	L
<b>Article Nr.</b>	3050-1	3050-2	3050-3



**CODE : 3025**

**FROG SPLINT**

*Used for finger fractures and deformities, sports injuries, and detect finger, and when parts of the finger need rest.*



<b>Finger Circ. (cm)</b>	4,5-5	5-6	6-7
<b>Size</b>	S	M	L
<b>Article Nr.</b>	3025-1	3025-2	3025-3



**CODE : 6500**

**KNEE PAD WITH ADJUSTABLE ANGLE**

*Used for internal, external and cross-link lesions the patella and femur traumas, also at knee joint when necessary to reduce to pressure. Produced from antiallergic and antibacterial sponge fabric.*



<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>6500</b>



**CODE : 7000**

**KNEE IMMOBILIZER**

*Provides support and stabilisation for leg and knee. Used for leg and knee for keep them stable when necessary. Produced of antiallergic and antibacterial sponge fabric.*



<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>7000</b>



**CODE : 6510**

**ROM BRACE**

Produced anatomically adjustable knee orthosis. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°-0°, 15°, 30°, 45°, 60°) angle adjustment possible. Immobilization limitation at (0°, 15°, 30°, 45°, 60°). Provides support and stability for knee injury. Used for protection at activity during post-op rehab following ligamentous surgical repairs and reconstruction. Controlled R.O.M. for post-traumatic, post surgical repair of human knee. It has one sizes. Length adjustment is possible from top and bottom.



Size	Standard
Article Nr.	6510



**CODE : 6520**

**ANGLE ADJUSTABLE KNEE PAD**

Produced anatomically adjustable knee orthosis. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°-0°, 15°, 30°, 45°, 60°) angle adjustment possible. Immobilization limitation at (0°, 15°, 30°, 45°, 60°) Provides support and stability for knee injury. Used for protection at activity during post-op rehab following ligamentous surgical repairs and reconstruction. Controlled R.O.M. for post-traumatic, post-surgical repair of human knee. It has one sizes. Length adjustment is possible from top and bottom.



Size	Standard
Article Nr.	6520



**CODE : 6530**

**KNEE PAD WITH ADJUSTABLE ANGLE**

Produced of light aluminum, functional, angle adjusted hinged rigid brace offers limitation on extension and flexion, shape of the knee. Flexion and extension angle adjustment is possible. Provides the most convenient stabilization with the 4-point principle and prevents tibial displacement and knee hyperextension.



Size	Standard
Article Nr.	6530



**CODE : 6540**

**KNEE PAD WITH ANGLE ADJUSTABLE HINGE**

Used for medial and lateral ligament instabilities, knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Angle adjustable hinge provides stabilisation in needed angle. Used for protection in sport activities. Produced from antiallergic, and antibacterial neoprene fabric.



Knee Circ. (cm)	30-34	34-38	38-42	42-46	46-50
Size	S	M	L	XL	XXL
Article Nr.	6540-S	6540-M	6540-L	6540-XL	6540-XXL



**CODE : 6550**

**KNEE PAD WITH ANGLE ADJUSTABLE HINGE ADVANCED**

Used for medial and lateral ligament instabilities, knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Angle adjustable hinge provides stabilisation in needed angle. Used for protection in sport activities. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6550-S</b>	<b>6550-M</b>	<b>6550-L</b>	<b>6550-XL</b>	<b>6550-XXL</b>



**CODE : 6302**

**LIGAMENT KNEE PAD WITH PATELLA TENDON SUPPORT**

Used for patella disorders, at arthritis, minor ligament insufficiency, ligament looseness, meniscus degenerations, and after operation. Support to lateral, and medial ligament with stays on both sides of knee pad. Offers protection with silicone pad in patella area. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6302-S</b>	<b>6302-M</b>	<b>6302-L</b>	<b>6302-XL</b>	<b>6302-XXL</b>



**CODE : 6402**

**HINGED KNEE PAD WITH PATELLA LIGAMENT SUPPORT**

Used for medial and lateral ligament instabilities, knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Offers protection with silicone pad in patella area. Used for protection in sport activities. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6402-S</b>	<b>6402-M</b>	<b>6402-L</b>	<b>6402-XL</b>	<b>6402-XXL</b>



**CODE : NV.6315**

**KNITTED KNEE PAD WITH HINGE**

Used for medial and lateral ligament instabilities, knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Offers protection with silicone pad in patella area. Used for protection in sport activities. Produced from antiallergic, and antibacterial knitted fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>NV.6315-S</b>	<b>NV.6315-M</b>	<b>NV.6315-L</b>	<b>NV.6315-XL</b>	<b>NV.6315-XXL</b>



**CODE : 6300**

**KNEE PAD WITH LIGAMENT SUPPORT**

Used for patella disorders, at arthritis, minor ligament insufficiency, ligament looseness, meniscus degenerations, and after operation. Support to lateral, and medial ligament with stays on both sides of knee pad. Produced from antiallergic, and antibacterial neoprene fabric.



Knee Circ. (cm)	30-34	34-38	38-42	42-46	46-50
Size	S	M	L	XL	XXL
Article Nr.	6300-S	6300-M	6300-L	6300-XL	6300-XXL



**CODE : 6400**

**KNEE PAD WITH HINGE**

Used for medial and lateral ligament instabilities, knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities. Produced from antiallergic, and antibacterial neoprene fabric.



Knee Circ. (cm)	30-34	34-38	38-42	42-46	46-50
Size	S	M	L	XL	XXL
Article Nr.	6400-S	6400-M	6400-L	6400-XL	6400-XXL



**CODE : 6200**

**KNEE PAD WITH PATELLA OPEN**

Used for knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6200-S</b>	<b>6200-M</b>	<b>6200-L</b>	<b>6200-XL</b>	<b>6200-XXL</b>




**CODE : 6250**

**KNEE PAD WITH VENCRO**

Used for knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities. Adjustable straps on top and bottom provides pressure over knee, and offers comfortable usage with providing maximum support to the knee.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6250-S</b>	<b>6250-M</b>	<b>6250-L</b>	<b>6250-XL</b>	<b>6250-XXL</b>



**CODE : 6150**

**PATELLA CLOSED KNEE PAD**

Used for knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6150-S</b>	<b>6150-M</b>	<b>6150-L</b>	<b>6150-XL</b>	<b>6150-XXL</b>



**CODE : 6301**

**KNEE PAD WITH ANTERIOR CRUCIATE SUPPORT**

Offers support with pad, and patella area is open to provide good fit. Provides extra stabilisation with spiral stays in joints. Offers stabilisation with 3 point cross straps to anterior cruciate ligament. Produced from antibacterial and antiallergic neoprene fabric.



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6301-S</b>	<b>6301-M</b>	<b>6301-L</b>	<b>6301-XL</b>	<b>6301-XXL</b>



**CODE : NV.6320**

**KNITTED KNEE PAD WITH LIGAMENT SUPPORT**

Used for knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities. Offers protection in patella area with silicone pad and offers stabilisation with stays. Produced from breathable, and flexible knitted fabric.



<b>Knee Circ. (cm)</b>	30-34	34-38	38-42	42-46	46-50
<b>Size</b>	S	M	L	XL	XXL
<b>Article Nr.</b>	NV.6320-S	NV.6320-M	NV.6320-L	NV.6320-XL	NV.6320-XXL



**CODE : NV.6323**

**KNITTED KNEE PAD CLASSIC**

Used for knee traumas, all inflammatory disorders in the knee joint, knee osteoarthritis, and stabilisation of patellar after operations. Used for protection in sport activities. Produced from breathable, and flexible knitted fabric.



<b>Knee Circ. (cm)</b>	30-34	34-38	38-42	42-46	46-50
<b>Size</b>	S	M	L	XL	XXL
<b>Article Nr.</b>	NV.6323-S	NV.6323-M	NV.6323-L	NV.6323-XL	NV.6323-XXL



**CODE : 800**

**WOOL KNEE PAD**

*Used for all kind of knee pain, rheumatism complaints, and keeping the knee warm. Produced anatomically, elastic wool fabric.*



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>800-S</b>	<b>800-M</b>	<b>800-L</b>	<b>800-XL</b>	<b>800-XXL</b>



**CODE : 810**

**KNEE PAD BASIC**

*The design is anatomically shaped to the knee. It has medium pressure. It is made from flexible neoprene fabric. It is suitable for daily activities and protect the knee from minor injuries.*



<b>Knee Circ. (cm)</b>	<b>30-34</b>	<b>34-38</b>	<b>38-42</b>	<b>42-46</b>	<b>46-50</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>810-S</b>	<b>810-M</b>	<b>810-L</b>	<b>810-XL</b>	<b>810-XXL</b>



**CODE : 6210**

**PATELLA TENDON BANDAGE**

*Anatomically designed, and produced from neoprene fabric. Has specially designed silicone pad that supports patella tendon. Does not prevent the movements of joint. Relieves tension and pain of tendon with its compression effect.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>6210</b>



**CODE : 9050**

**ROM WALKER**

*Used for foot and ankle traumas, tendon ligament soft tissue injuries and fractures also it is for stabilisation after operations.*



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9050-S	9050-M	9050-L



**CODE : 9060**

**ROM WALKER (SHORT)**

*Used for foot and ankle traumas, tendon ligament soft tissue injuries and fractures also it is for stabilisation after operations.*



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9060-S	9060-M	9060-L



**CODE : 9070**

**ACHILL BOOT**

Used for correction and stabilisation of soft tissue for the large toe. Internal pad produced from terry cloth. Offers stabilisation in needed position. Gradual setting on the sole offers different position for foot.



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9070-S	9070-M	9070-L



**CODE : 9080**

**ACHILL BOOT (SHORT)**

Used for correction and stabilisation of soft tissue for the large toe. Internal pad produced from terry cloth. Offers stabilisation in needed position. Gradual setting on the sole offers different position for foot.



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9080-S	9080-M	9080-L



**CODE : NV.3310**

**KNITTED ANKLE PAD WITH SILICONE SUPPORT**

Used for ankle sprains to stabilization of ligament, soft tissue injuries. Suitable for all sport activities. Provide support to ankle with silicone pads. Produced from breathable, and flexible knitted fabric.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	NV.3310-S	NV.3310-M	NV.3310-L	NV.3310-XL	NV.3310-XXL



**CODE : NV.3313**

**KNITTED ANKLE PAD WITH SILICONE SUPPORT WITH STRAPS**

Used for ankle sprains to stabilization of ligament, soft tissue injuries. Suitable for all sport activities. Provide support to ankle with silicone pads. Offers individual stabilisation with straps. Produced from breathable, and flexible knitted fabric.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	NV.3313-S	NV.3313-M	NV.3313-L	NV.3313-XL	NV.3313-XXL



**CODE : NV.3320**

**KNITTED ANKLE PAD WITH ACHILL TENDON SUPPORT**

Used for ankle sprains to stabilization of ligament, soft tissue injuries. Suitable for all sport activities. Provide support to ankle with silicone pads in under heel and on ankle. Produced from breathable, and flexible knitted fabric.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	NV.3320-S	NV.3320-M	NV.3320-L	NV.3320-XL	NV.3320-XXL



**CODE : NV.3305**

**KNITTED ANKLE PAD CLASSIC**

Used for ankle traumas to relieve the pain, and to support the ankle. Preserve the body temperature due to neoprene material. Accelerates the recovery process by improving blood circulation. Produced from antiallergic, antibacterial knitted fabric.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	NV.3305-S	NV.3305-M	NV.3305-L	NV.3305-XL	NV.3305-XXL



**CODE : 6810**

**ANKLE PAD SIMPLE**

*Used for ankle traumas to relieve the pain, and to support the ankle. Preserve the body temperature due to neoprene material. Accelerates the recovery process by improving blood circulation. Produced from antiallergic, antibacterial neoprene fabric.*



<b>Ankle Circ. (cm)</b>	<b>21-23</b>	<b>23-25</b>	<b>25-27</b>	<b>27-29</b>	<b>29-31</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6810-S</b>	<b>6810-M</b>	<b>6810-L</b>	<b>6810-XL</b>	<b>6810-XXL</b>



**CODE : 9040**

**DORSIFLEXION BANDAGE**

*Produced from sponge fabric. Offers needed tension for foot and use for low foot deformation dorsiflexion. Provides posture correction of foot for patients. Available for usage inside of shoes.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>9040</b>



**CODE : 3300**

**ANKLE PAD WITH SILICONE SUPPORT**

Used in ankle traumas to support the ankle and relieve the pain. Offers support with silicone pad. Used for supporting ankle sprains and injuries and to support the ankle. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Ankle Circ. (cm)</b>	<b>21-23</b>	<b>23-25</b>	<b>25-27</b>	<b>27-29</b>	<b>29-31</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>3300-S</b>	<b>3300-M</b>	<b>3300-L</b>	<b>3300-XL</b>	<b>3300-XXL</b>



**CODE : 9000**

**ANKLE PAD WITH LIGAMENT SUPPORT**

Used in ankle traumas to support the ankle and relieve the pain. Has 8 shaped of straps which wraps the ankle. Used for supporting ankle sprains and injuries and to support the ankle. Produced from antiallergic, and antibacterial neoprene fabric.



<b>Ankle Circ. (cm)</b>	<b>21-23</b>	<b>23-25</b>	<b>25-27</b>	<b>27-29</b>	<b>29-31</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>9000-S</b>	<b>9000-M</b>	<b>9000-L</b>	<b>9000-XL</b>	<b>9000-XXL</b>



**CODE : 9025**

**REFLEX AFO**

Used for fixing ankle at 90 degree angle also used for nerve and foot traumas. Right and left sizes options are available.



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9025-1 9025-S1	9025-2 9025-S2	9025-3 9025-S3



**CODE : 9030**

**FOOT GUARD - T KIT**

Produced from thermoplastic material. Parts where touches the skin is covered with terry cloth. Offers stabilisation of foot at 90°. Designed as T kit which prevents rotation of foot.



Foot Size	35-38	39-42	43-46
Size	S	M	L
Article Nr.	9030-S	9030-M	9030-L



**CODE : 9150**

### **AIR CAST WITH PAD**

*Produced from special thermo plastic material. Offers relaxation with pad. Used for fixing ankle. Easy to use with maximum ankle support.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>9150</b>



**CODE : 9160**

### **AIR CAST**

*Produced from special thermo plastic material. Offers relaxation with air pad. Used for fixing ankle. Easy to use with maximum ankle support.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>9160</b>



**CODE : 9170**

### **AIR CAST WITH GEL**

*Produced from special thermo plastic material. Offers relaxation with gel pad. Used for fixing ankle. Easy to use with maximum ankle support.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>9170</b>



**CODE : 9001**

**ANKLE PAD WITH LIGAMENT SUPPORT ELEGANCE**

Used for ankle traumas, ankle sprains, injuries and relieves pain, and aiding for recovery process.

Produced from antiallergic, and antibacterial neoprene fabric.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	9001-S	9001-M	9001-L	9001-XL	9001-XXL



**CODE : 2900**

**METATARSAL PAD**

Used for when metatars curvature spoiled and aids to reform the metatars curvature.

Size	Standard
Article Nr.	2900



**CODE : 5000**

**HALLUX VALGUS NIGHT SPLINT**

*Used to keep thumb bent outwards at hallux valgus problems. Used for treatment after operation for recovery process. Suitable for using at nights. Produced from sponge fabric.*

<b>Size</b>	<b>Right</b>	<b>Left</b>
<b>Article Nr.</b>	<b>5000-R</b>	<b>5000-L</b>



**CODE : 6800**

**ANKLE PAD WITH VENCRO**

*Used for ankle traumas, and for relieving pain. Produced from antiallergic and antibacterial neoprene fabric.*



<b>Ankle Circ. (cm)</b>	<b>21-23</b>	<b>23-25</b>	<b>25-27</b>	<b>27-29</b>	<b>29-31</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>6800-S</b>	<b>6800-M</b>	<b>6800-L</b>	<b>6800-XL</b>	<b>6800-XXL</b>



**CODE : 5210**

**CORK INSOLE**

*Provides support and reduce stress in the joints by spreading body burden equally when standing. Offers longitudinal and metatarsal support. Produced from special thermoplastic material.*

Foot Size	33	34	35	36	37	38	39	40	41	42	43	44	45	46
Article Nr.	5210-33 to 46													




**CODE : 5200**

**STEEL INSOLE**

*Recommended to use in the treatment of pes cavus, pes planus, epin calcainen and metatarsal deformations. Produced from leather and supported with steel bridge.*

Foot Size	33	34	35	36	37	38	39	40	41	42	43	44	45	46
Article Nr.	5200-33 to 46													



**CODE : 7400**

**UPPER CALF PAD**

*Used for fiber breaks and injuries, spasms, in physical therapy process. Used for protective in sport activities. Produced from antibacterial and antiallergic neoprene fabric.*



Upper Calf Circ. (cm)	45-48	48-52	52-56	56-60	60-64
Size	S	M	L	XL	XXL
Article Nr.	7400-S	7400-M	7400-L	7400-XL	7400-XXL



**CODE : NV.7401**

**KNITTED UPPER CALF PAD**

*Used for fiber breaks and injuries, spasms, in physical therapy process. Used for protective in sport activities. Produced from breathable, flexible, antibacterial and antiallergic knitted fabric.*



Upper Calf Circ. (cm)	45-48	48-52	52-56	56-60	60-64
Size	S	M	L	XL	XXL
Article Nr.	NV.7401-S	NV.7401-M	NV.7401-L	NV.7401-XL	NV.7401-XXL




**CODE : 7500**

**LOWER CALF PAD**

Used for lower fiber breaks and injuries, spasms, in process of physical therapy and also various sportactivities as a protector. Produced from antiallergic and antibacterial neoprene fabric.



Lower Calf Circ. (cm)	27-30	30-33	33-36	36-39	39-42
Size	S	M	L	XL	XXL
Article Nr.	7500-S	7500-M	7500-L	7500-XL	7500-XXL



**CODE : NV.7501**

**KNITTED LOWER CALF PAD**

Used for lower fiber breaks and injuries, spasms, in process of physical therapy and also various sportactivities as a protector. Produced from breathable, flexible, antiallergic and antibacterial knitted fabric.



Lower Calf Circ. (cm)	27-30	30-33	33-36	36-39	39-42
Size	S	M	L	XL	XXL
Article Nr.	NV.7501-S	NV.7501-M	NV.7501-L	NV.7501-XL	NV.7501-XXL



**CODE : 7900**

**THERMAL SHORTS**

*Increase blood circulation, and keeps body heat. Prevent occurring edema with applied compression. Relieves pain by reducing tension in muscles. Designed for protection to inguen damages. Produced from antiallergic, antibacterial neoprene fabric.*



<b>Waist Circ. (cm)</b>	<b>70-80</b>	<b>80-90</b>	<b>90-100</b>	<b>100-110</b>	<b>110-120</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>7900-S</b>	<b>7900-M</b>	<b>7900-L</b>	<b>7900-XL</b>	<b>7900-XXL</b>



**CODE : 7901**

**HIP SUPPORT**

*Produced from antibacterial and antiallergic neoprene fabric. Supports blood circulation, and keeps body heat. Designed as specially to anatomical shape of hip. Provides compression to the upper calf and hip joint by shaped support. Offers protection and reduce muscle tensions in hip and upper calf area.*



<b>Upper Calf Circ. (cm)</b>	<b>45-48</b>	<b>48-52</b>	<b>52-56</b>	<b>56-60</b>	<b>60-64</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>	<b>XXL</b>
<b>Article Nr.</b>	<b>7901-1</b> <b>7901-S1</b>	<b>7901-2</b> <b>7901-S2</b>	<b>7901-3</b> <b>7901-S3</b>	<b>7901-4</b> <b>7901-S4</b>	<b>7901-5</b> <b>7901-S5</b>



**CODE: NV.6100**

**ELASTIC KNEE PAD**

Used for muscle strains, and after rehabilitation after injuries during sport activities.



Knee Circ. (cm)	30-34	34-38	38-42	42-46	46-50
Size	S	M	L	XL	XXL
Article Nr.	NV.6100-S	NV.6100-M	NV.6100-L	NV.6100-XL	NV.6100-XXL



**CODE: 7600**

**ELASTIC LOWER CALF PAD**

Used for supporting to front side of leg with applied compression. Used for muscle strains, and injuries during sport activities.



Lower Calf Circ.(cm)	27-30	30-33	33-36	36-39	39-42
Size	S	M	L	XL	XXL
Article Nr.	7600-S	7600-M	7600-L	7600-XL	7600-XXL



**CODE: NV.6321**

**ELASTIC ANKLE PAD**

Ideal for enforced and fatigued ankles. Provides compression to necessary areas and allows for motions. Especially used in sport activities.



Ankle Circ. (cm)	21-23	23-25	25-27	27-29	29-31
Size	S	M	L	XL	XXL
Article Nr.	NV.6321-S	NV.6321-M	NV.6321-L	NV.6321-XL	NV.6321-XXL



## **PILLOWS**

*Orthopedic Visco pillow of Novamed are necessary to use for neck hernia, at arm numbness, patients with neurological problems and that also used for snoring problems which is caused by incorrect sleeping position. Because of anatomical shape, it avoid possible causes of neck and shoulder pains. Produced from visco memory foam.*

**VISCO PILLOW**



**Small (50\*30)      Large (60\*40)**  
**Code : 120-K      Code : 120-B**

**MEDICAL VISCO PILLOW**



**Code : SR-1280**

**FAT VISCO PILLOW**



**Code : SR-1281**

**REFLUX PILLOW**



**Code : SR-1191**

**TRAVEL PILLOW**



**Code : SR-1192**

**SNORING PILLOW**



**Code : SR-1193**

**WAIST SUPPORT**

*Code : 140*


**SITTING PAD**

*Code : 130*

**SEATING CUSHION**

*Code : SR-1279*

**Medical Stocking Size Table**



Size (cm)							
Leg Circ.	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	22-26	24-29	27-33	29-36	31-38	33-39	35-40
cC	28-33	30-37	33-40	35-43	37-45	39-46	41-48
cD	27-32	28-35	32-39	34-42	36-44	38-45	40-47
cE	32-38	33-41	37-44	39-47	41-50	43-52	45-54
cF	37-47	38-52	43-58	46-61	50-66	53-69	56-71
cG	43-57	46-60	49-67	52-72	55-76	60-79	64-80

Length of leg (cm)	
<u>lG</u>	71-82
<u>lG</u>	62-70
<u>lD</u>	38-43
<u>lD</u>	33-37

**Lymph Edema Arm Sleeve Size Table**





Size (cm)							
Arm Circ.	1	2	3	4	5	6	7
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cE	16-22	21-27	24-30	27-33	30-36	33-39	36-42
cG	20-27	23-30	26-33	30-37	34-41	38-45	42-49


Length of arm (cm)	
<u>lG</u>	Over 44 cm
<u>lG</u>	Up to 44 cm




**CCL1 UNDER KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 18-21mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.600  
Closed Toe

 Code : NV.610  
Open Toe


Size	1	2	3	4	5	6	7
Article Nr.	NV.600-1	NV.600-2	NV.600-3	NV.600-4	NV.600-5	NV.600-6	NV.600-7
	NV.610-1	NV.610-2	NV.610-3	NV.610-4	NV.610-5	NV.610-6	NV.610-7




**CCL2 UNDER KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 23-32mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.660  
Closed Toe



 Code : NV.670  
Open Toe


Size	1	2	3	4	5	6	7
Article Nr.	NV.660-1	NV.660-2	NV.660-3	NV.660-4	NV.660-5	NV.660-6	NV.660-7
	NV.670-1	NV.670-2	NV.670-3	NV.670-4	NV.670-5	NV.670-6	NV.670-7



**CCL3 UNDER KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 34-46mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.720  
Open Toe


Size	1	2	3	4	5	6	7
Article Nr.	NV.720-1	NV.720-2	NV.720-3	NV.720-4	NV.720-5	NV.720-6	NV.720-7




**CCL1 ABOVE KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It avoids slipping by silicone band. It treats edema and varicosis by 18-21mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option;it gives perfect result for the diabetic patients,patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.620  
Closed Toe

 Code : NV.630  
Open Toe


Size	1	2	3	4	5	6	7
Article Nr.	NV.620-1	NV.620-2	NV.620-3	NV.620-4	NV.620-5	NV.620-6	NV.620-7
	NV.630-1	NV.630-2	NV.630-3	NV.630-4	NV.630-5	NV.630-6	NV.630-7




**CCL2 ABOVE KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It avoids slipping by silicone band. It treats edema and varicosis by 23-32mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option;it gives perfect result for the diabetic patients,patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.680  
Closed Toe

 Code : NV.690  
Open Toe


Size	1	2	3	4	5	6	7
Article Nr.	NV.680-1	NV.680-2	NV.680-3	NV.680-4	NV.680-5	NV.680-6	NV.680-7
	NV.690-1	NV.690-2	NV.690-3	NV.690-4	NV.690-5	NV.690-6	NV.690-7



**CCL3 ABOVE KNEE COMPRESSION STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It avoids slipping by silicone band.It treats edema and varicosis by 34-46mmHg pressure also prevents swellings at leg. It is able to use for beginning of varicosis and after the leg, vein operations. By the open toe option;it gives perfect result for the diabetic patients,patient with nutritional deficiency at skin and patient with toe deformation.

Colors  

 Code : NV.730  
Open Toe

Size	1	2	3	4	5	6	7
Article Nr.	NV.730-1	NV.730-2	NV.730-3	NV.730-4	NV.730-5	NV.730-6	NV.730-7



### CCL1 PANTYHOSE COMPRESSION STOCKING

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 18-21mmHg pressure also prevents swellings at leg. It is ideal for usage for beginning of varicosis and after vein operation at the legs. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin, and patient with toe deformation.

Colors



Code : NV.640  
Closed Toe



Code : NV.650  
Open Toe

Size	1	2	3	4	5	6	7
Article Nr.	NV.640-1 NV.650-1	NV.640-2 NV.650-2	NV.640-3 NV.650-3	NV.640-4 NV.650-4	NV.640-5 NV.650-5	NV.640-6 NV.650-6	NV.640-7 NV.650-7



### CCL2 PANTYHOSE COMPRESSION STOCKING

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 23-32mmHg pressure also prevents swellings at leg. It is ideal for usage for beginning of varicosis and after vein operation at the legs. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin, and patient with toe deformation.

Colors



Code : NV.700  
Closed Toe



Code : NV.710  
Open Toe

Size	1	2	3	4	5	6	7
Article Nr.	NV.700-1 NV.710-1	NV.700-2 NV.710-2	NV.700-3 NV.710-3	NV.700-4 NV.710-4	NV.700-5 NV.710-5	NV.700-6 NV.710-6	NV.700-7 NV.710-7



### CCL3 PANTYHOSE COMPRESSION STOCKING

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis by 34-46mmHg pressure also prevents swellings at leg. It is ideal for usage for beginning of varicosis and after vein operation at the legs. By the open toe option; it gives perfect result for the diabetic patients, patient with nutritional deficiency at skin, and patient with toe deformation.

Colors




Code : NV.740  
Open Toe

Size	1	2	3	4	5	6	7
Article Nr.	NV.740-1	NV.740-2	NV.740-3	NV.740-4	NV.740-5	NV.740-6	NV.740-7



**UNDER KNEE ANTI EMBOLISM STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It is used before or after operations also used for cellulite, gangrene dermatitis, edema and skin infections.


 Code : NV.755  
Open Toe

Size	1	2	3	4	5
Article Nr.	NV.755-1	NV.755-2	NV.755-3	NV.755-4	NV.755-5



**ABOVE KNEE ANTI EMBOLISM STOCKING**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It avoids slipping by silicone band. It is used before or after operations also used for cellulite, gangrene dermatitis, edema and skin infections.


 Code : NV.750  
Open Toe

Size	1	2	3	4	5
Article Nr.	NV.750-1	NV.750-2	NV.750-3	NV.750-4	NV.750-5



**ANTI EMBOLISM STOCKING WITH WAIST BELT**

It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It is used before or after operations also used for cellulite, gangrene dermatitis, edema and skin infections.

 Code : NV.745  
Open Toe

Size	1	2	3	4	5
Article Nr.	NV.745-1	NV.745-2	NV.745-3	NV.745-4	NV.745-5



**PREGNANT COMPRESSION STOCKING**

*It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It treats edema and varicosis also prevents swellings at leg. It is able to use for during pregnancy and after given birth .*



**Code : NVH.500**  
**Closed Toe**

<b>Size</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>
<b>Article Nr.</b>	<b>NVH.500-1</b>	<b>NVH.500-2</b>	<b>NVH.500-3</b>	<b>NVH.500-4</b>	<b>NVH.500-5</b>	<b>NVH.500-6</b>	<b>NVH.500-7</b>



**TRAVEL STOCKING**

*It is knitted with 1st class yarn by newest technology. It provides comfortable use with elastic knitting technic and gives all day rest to the legs. It is used for long journeys. Accelerates blood circulation for immobile legs. It prevents to swellings at legs also provides healthy and comfortable journey.*



**Code : NV.770**  
**Closed Toe**

<b>Size</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
<b>Article Nr.</b>	<b>NV.770-1</b>	<b>NV.770-2</b>	<b>NV.770-3</b>	<b>NV.770-4</b>	<b>NV.770-5</b>



**CCL1 LYMPH EDEMA ARM SLEEVE (CLASSIC)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



**Code : NV.760**  
**With Finger**

**Code : NV.761**  
**Without Finger**

Size	1	2	3	4	5	6	7
Article Nr.	NV.760-1	NV.760-2	NV.760-3	NV.760-4	NV.760-5	NV.760-6	NV.760-7
	NV.761-1	NV.761-2	NV.761-3	NV.761-4	NV.761-5	NV.761-6	NV.761-7



**CCL2 LYMPH EDEMA ARM SLEEVE (CLASSIC)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



**Code : NV.762**  
**With Finger**

**Code : NV.763**  
**Without Finger**

Size	1	2	3	4	5	6	7
Article Nr.	NV.762-1	NV.762-2	NV.762-3	NV.762-4	NV.762-5	NV.762-6	NV.762-7
	NV.763-1	NV.763-2	NV.763-3	NV.763-4	NV.763-5	NV.763-6	NV.763-7



**CCL3 LYMPH EDEMA ARM SLEEVE (CLASSIC)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



**Code : NV.764**  
**With Finger**

**Code : NV.765**  
**Without Finger**

Size	1	2	3	4	5	6	7
Article Nr.	NV.764-1	NV.764-2	NV.764-3	NV.764-4	NV.764-5	NV.764-6	NV.764-7
	NV.765-1	NV.765-2	NV.765-3	NV.765-4	NV.765-5	NV.765-6	NV.765-7



**CCL1 LYMPH EDEMA ARM SLEEVE (SILICONE BAND)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. It also avoids slipping by silicone band. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.

Colors  

Code : NV.781  
With Finger

Code : NV.766  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.781-1	NV.781-2	NV.781-3	NV.781-4	NV.781-5	NV.781-6	NV.781-7
	NV.766-1	NV.766-2	NV.766-3	NV.766-4	NV.766-5	NV.766-6	NV.766-7



**CCL2 LYMPH EDEMA ARM SLEEVE (SILICONE BAND)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. It also avoids slipping by silicone band. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.

Colors  

Code : NV.767  
With Finger

Code : NV.768  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.767-1	NV.767-2	NV.767-3	NV.767-4	NV.767-5	NV.767-6	NV.767-7
	NV.768-1	NV.768-2	NV.768-3	NV.768-4	NV.768-5	NV.768-6	NV.768-7



**CCL3 LYMPH EDEMA ARM SLEEVE (SILICONE BAND)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. It also avoids slipping by silicone band. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.

Colors  

Code : NV.769  
With Finger

Code : NV.779  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.769-1	NV.769-2	NV.769-3	NV.769-4	NV.769-5	NV.769-6	NV.769-7
	NV.779-1	NV.779-2	NV.779-3	NV.779-4	NV.779-5	NV.779-6	NV.779-7



**CCL1 LYMPH EDEMA ARM SLEEVE (WITH STRAP)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. Shoulder strap support is available. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



Code : NV.771  
With Finger

Code : NV.772  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.771-1 NV.772-1	NV.771-2 NV.772-2	NV.771-3 NV.772-3	NV.771-4 NV.772-4	NV.771-5 NV.772-5	NV.771-6 NV.772-6	NV.771-7 NV.772-7



**CCL2 LYMPH EDEMA ARM SLEEVE (WITH STRAP)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. Shoulder strap support is available. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



Code : NV.773  
With Finger

Code : NV.774  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.773-1 NV.774-1	NV.773-2 NV.774-2	NV.773-3 NV.774-3	NV.773-4 NV.774-4	NV.773-5 NV.774-5	NV.773-6 NV.774-6	NV.773-7 NV.774-7



**CCL3 LYMPH EDEMA ARM SLEEVE (WITH STRAP)**

It allows to make condition of damaged lymphatic vein with its compression and it avoids stop up blood circulation after chest surgeries. In production high quality and special materials are used. Shoulder strap support is available. These materials are anti allergic and environmentally harmless. Arm sleeves with finger and without finger is available for all types.



Code : NV.775  
With Finger

Code : NV.778  
Without Finger

Size	1	2	3	4	5	6	7
Article Nr.	NV.775-1 NV.778-1	NV.775-2 NV.778-2	NV.775-3 NV.778-3	NV.775-4 NV.778-4	NV.775-5 NV.778-5	NV.775-6 NV.778-6	NV.775-7 NV.778-7

## **SILICONE PRODUCTS**

*Novamed Silicone Products with innovative models cover your needs, not just in sport activities but also in daily life.*

*Latest models of silicone products with gel technology, Novamed offers best products.*



**CODE : 100**

**SILICONE INSOLE**

*Provides balanced distribution of body burden of standing, and reduces the burden and stress of joints. Used at standing for long time. Produced from silicone in accordance with anatomical shape of the foot.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	100-0	100-1	100-2	100-3	100-4	100-5



**CODE : 125**

**SILICONE INSOLE COVERED WITH FABRIC**

*Provides balanced distribution of body burden of standing, and reduces the burden and stress of joints. Used at standing for long time. Recommended for employees who work as standing. Produced from silicone in accordance with anatomical shape of the foot.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	125-0	125-1	125-2	125-3	125-4	125-5



**CODE : 402**

**SILICONE PERFORATED INSOLE**

*Enables even distribution of pressure on the sole by providing anatomical support. Has relaxing and comforting effect.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	402-0	402-1	402-2	402-3	402-4	402-5



**CODE : 101**

**SILICONE PERFORATED DOTTED INSOLE**

*Offers even distribution of pressure on the sole by providing anatomical support. Presents relaxing and comforting effect. Offers an anti slip feature by adhering to the foot.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	101-0	101-1	101-2	101-3	101-4	101-5



**CODE : 181**

**SILICONE INSOLE ARC SUPPORT**

*Provides balanced distribution of body burden of standing, and reduces the burden and stress of joints. Used at standing for long time. Recommended for employees who work as standing. Produced from silicone in accordance with anatomical shape of the foot.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	181-0	181-1	181-2	181-3	181-4	181-5



**CODE : 115**

**SILICONE HALF INSOLE**

*Maximum elastic component, hypoallergenic, and orthopedic. Used in heel spur treatment. Prevent waist and leg pain which occur for standing long time, and foot pain. Recommended for employees who work as standing. Produced from silicone in accordance with anatomical shape of the foot.*



Foot No.	35-36	37-38	39-40	41-42	43-44	45-46
Size	XS	S	M	L	XL	XXL
Article Nr.	115-0	115-1	115-2	115-3	115-4	115-5



**CODE : 110**

**SILICONE HEEL**

*Used in heel spur treatment. Reduces pain over heel which is caused by burden. Produced from silicone in accordance with anatomical shape of the foot.*



<b>Foot Size</b>	<b>35-38</b>	<b>39-42</b>	<b>43-46</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>110-1</b>	<b>110-2</b>	<b>110-3</b>



**CODE : 111**

**SILICONE KAILE 3 MM HEEP CUP**

*Compensates for difference in case of length gap between legs.*



<b>Foot No.</b>	<b>35-37</b>	<b>38-40</b>	<b>41-43</b>	<b>44-46</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
<b>Article Nr.</b>	<b>111-1</b>	<b>111-2</b>	<b>111-3</b>	<b>111-4</b>



**CODE : 4220**

**SILICONE ADHESIVE HEEL CUP**

*Used for plantar fasciitis, calcaneal pains and heel spur treatments. Available for usage in all kinds of shoes.*



Foot No.	35-37	38-40	41-43	44-46
Size	S	M	L	XL
Article Nr.	4220-1	4220-2	4220-3	4220-4



**CODE : 4240**

**SILICONE ADHESIVE METATARSAL PAD**

*Provides relieve pain in metatarsal area with pressure.*



Foot No.	35-40	41-46
Size	S	M
Article Nr.	4240-1	4240-2



**CODE : 180**

**SILICONE ARC SUPPORT**

*Provides advanced support for arches with effective gel based. Offers maximum relieve from foot strain. Used for flat feed, rigid steps, and high or weak fallen arches.*



<b>Foot Size</b>	<b>35-36</b>	<b>37-38</b>	<b>39-40</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>180-1</b>	<b>180-2</b>	<b>180-3</b>



**CODE : 5215**

**SILICONE MEDIAL AND LATERAL WEDGE**

*Produced from silicone material, and covered with fabric. Placed on inner and outside of the shoe. Suitable for daily use in shoes.*



<b>Foot No.</b>	<b>35-37</b>	<b>38-40</b>	<b>41-43</b>	<b>44-46</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
<b>Article Nr.</b>	<b>5215-1</b>	<b>5215-2</b>	<b>5215-3</b>	<b>5215-4</b>



**CODE : 155**

**SILICONE TOE CREST WITH RING**

*Protects toes and prevent them from twisting. Helps to reduce pain and to treat toe ailments.*

<b>Size</b>	<b>Standard / Right - Left</b>
<b>Article Nr.</b>	<b>155-L 155-R</b>



**CODE : 417**

**SILICONE FOREFOOT INSOLE**

*Recommended as reliever to people that stands very long periods and especially to women wear high heels.*



<b>Size</b>	<b>Standard / Right - Left</b>
<b>Article Nr.</b>	<b>417-R 417-L</b>



**CODE : 5900**

**SILICONE PULLEY**

*Used for correction finger posture disorder with silicone attached between thumb and second finger. Suitable for daily use. Produced from high quality silicone.*



<b>Foot No.</b>	<b>35-38</b>	<b>38-42</b>	<b>43-46</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>5900-1</b>	<b>5900-2</b>	<b>5900-3</b>



**CODE : 165**

**SILICONE FINGER SEPERATOR**

*Used for reducing friction ay edges and overlapped fingers with silicone seperators. Reduces pain and irritation between fingers. Produced from high quality silicone.*



<b>Foot No.</b>	<b>35-38</b>	<b>38-42</b>	<b>43-46</b>
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>165-1</b>	<b>165-2</b>	<b>165-3</b>



**CODE : 175**

**SILICONE HEEL BANDAGE**

*Produced from hypoallergic sticker material. Silicone heel band includes %100 mineral oil. Reduces pain and rubbing.*

<b>Size</b>	<b>Standard</b>
-------------	-----------------

<b>Article Nr.</b>	<b>175</b>
--------------------	------------



**CODE : 160**

**SILICONE SHAFT METATARSAL PAD**

*Passes from the tip of the finger and placed to the metatarsal region. Used for metatarsalgia, hyperkeretosis and as a shock absorber. Produced from high quality silicone.*

<b>Size</b>	<b>Standard</b>
-------------	-----------------

<b>Article Nr.</b>	<b>160</b>
--------------------	------------



**CODE : 5220**

**SILICONE FINGER COVER**

*It is used for all kind of problems about ote fingers and to protect and relieve diabetes patients.*



Foot No.	35-38	38-42	43-46
Size	S	M	L
Article Nr.	5220-1	5220-2	5220-3



**CODE : 188**

**SILICONE TOE CAP**

*Helps proper positioning of overlapping big toe in the shoe.*



Foot No.	35-40	41-46
Size	S	L
Article Nr.	188-1	188-3



**CODE : 185**

**SILICONE BUNION PAD**

*Used in metatarsalgia and hyperkeratosis as absorvation to shock. Produced from high quality silicone.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>185</b>



**CODE : 186**

**SILICONE BUNION PAD WITH WEDGED**

*Used in metatarsalgia and hyperkeratosis as absorvation to shock. Also wedged silicone is available between toe and second finger. Produced from high quality silicone.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>186</b>



**CODE : 420**

**SILICONE METATARSAL PAD WITH SEPERATOR**

*Used as protector and supporter in metatarsal area. Relieves pain and prevents hardened areas from expanding. Designed for flip flops and open toe shoes.*

<b>Size</b>	<b>Standard / Right - Left</b>
<b>Article Nr.</b>	<b>420-R 420-L</b>



**CODE : 4130**

**SILICONE TO CREST WITH RING**

*Used for protection toes and prevent from twisting. Offers fast healing process and helps for reducing pain in toes.*

<b>Size</b>	<b>Standard / Right - Left</b>
<b>Article Nr.</b>	<b>4130-R 4130-L</b>



## **SUPPORT PRODUCTS**

*Novamed Support Products support and protect the musculoskeletal system. They contribute to the healing process in cases of musculoskeletal disorders and injuries. Novamed Products help people to increase their mobility, and used to stabilize, immobilize and relieve affected areas of the body. Supports to regain functions and prevent or reduce pain during persistent discomfort.*



**CODE : 3600**

**WALKER (FOLDING - FIX)**

*Produced from aluminium pipe. Moves back and forth while walking due to joint connection. It can be foldable and stored easily without spending space.*

**Size**                      **Standard**

**Article Nr.**            **3600**



**CODE : 3700**

**CLOSET EXTENDER (COMMUNE CHAIRS)**

*Produced from aluminium pipe and plastic where place of sitting. Can be placed onto toilet seat when potty plugged in the patient may resolve the need for toilet.*

**Size**                      **Standard**

**Article Nr.**            **3700**



**CODE : 3725**

**COMMODE CHAIRS**

*At the place of sitting has toilet seat cover, soft pad and the sliders. Also it can place onto toilet seat in order to use it more comfortable.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>3725</b>



**CODE : 3400**

**CRUTCH (DOUBLE)**

*It is made of aluminum tubing. Underarm and hand grip are produce of hard plastics.*

<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>
<b>Article Nr.</b>	<b>3400-1</b>	<b>3400-2</b>	<b>3400-3</b>



**CODE : 3500**

**FOREARM CRUTCH (DOUBLE)**

*Produced from metal material. Provides extra support for walk. Adjustable sizes are available.*

<b>Size</b>	<b>Standard</b>
-------------	-----------------

<b>Article Nr.</b>	<b>3500</b>
--------------------	-------------



**CODE : 3950**

**CANE**

*Produced from metal material. Provides extra support for walk. Adjustable sizes are available.*

<b>Size</b>	<b>Standard</b>
-------------	-----------------

<b>Article Nr.</b>	<b>3950</b>
--------------------	-------------



**CODE : 3875**

**QUAD BALANCED CANE**

*The bariatric quad canes are made out of aluminum and have a small base and an offset handle. They have a lightweight design and low center of gravity, which combine for greater stability and balance.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>3875</b>



**CODE : 3850**

**TRIO BALANCED CANE**

*It is made of metal and gives to extra support for walk. It produced one size and it is adjustable to all sizes.*

<b>Size</b>	<b>Standard</b>
<b>Article Nr.</b>	<b>3850</b>



**OTHER PRODUCTS**



**CODE : 4150**

**HOT COLD GEL**

*Offers hot and cold treatment in one product. Easy to use. Multi usable. Can be used both hot and cold. Offers to keep its temperature hot/cold for long time. Gel compression packages provides natural treatment in both minor and big injuries, and pain.*

Size	S (13*13cm)	M (13*26cm)	L (26*35cm)
Article Nr.	4150-1	4150-2	4150-3



**CODE : SB**

**ELASTIC BANDAGE**

Size	6cm	8cm	10cm	12cm	15cm	20cm
Article Nr.	SB-690	SB-695	SB-700	SB-705	SB-710	SB-715



**CODE : NAO**

**ALEZE PAD**

Size	90*190cm	160-190cm
Article Nr.	NAO-T	NAO-Ç





# NOVAMED®

*Your Health is important to us*

*Çınar Medikal Sağlık Hizmetleri  
Şanayi ve Ticaret Ltd.Şti.*

*Esenkent Mah. Enbiya Sok.  
No:25 Ümraniye / İSTANBUL  
Tel: +90 216 573 92 26  
Fax: +90 216 469 40 29*

*export@novamed.tc  
www.novamed.tc*